

Nursing Cart

Nursing carts are flexibly designed based on the nursing needs of medical institutions, combining functionality, specifications, and style. Clinical staff can precisely position nursing carts containing essential supplies and equipment.

Features



Compatible with ISO 6040 baskets/trays



Fine medicine management, orderly placement



The Lifting Overbridge can be equipped with accessories such as tissue box, glove box, and storage shelf.



· Medical sharps box



· Glove Box Holder

Item No: T003

Worktop: 651mm x 439mm x 46mm

Dimension: 754mm x 541mm x 1005mm

Package Material: Cardboard Box

N.W: 32.00KG