

Pediatric Emergency Crash Carts

The Broselow system is a color-coded system that assists in providing rapid and accurate dosages of medications for pediatric patients through color variations.

Features



Compatible with ISO 6040 baskets/trays





Fine medicine management, orderly placement



Utilizes the Broselow color-coding design system for quick dosage identification.



Equipped with rail locks to prevent drawers from sliding out during movement.



· CPR Board



· Medical sharps box

Item No: T007

Worktop: 651mm x 439mm x 46mm

Dimension: 754mm x 541mm x 1005mm

Package Material: Cardboard Box

N.W: 66.50KG